

## MOLECULAR BIOLOGY PROTOCOLS

### I. Tissue Homogenization & RNA Isolation for cDNA synthesis

1. Pre-spin 2ml Phase lock Gel-heavy tubes briefly (ie- 3,000rpm for 2 min) @ RT.
2. **HOMOGENIZATION:** Obtain two rat eyes from  $-80^{\circ}\text{C}$  freezer and place on ice in a clean uncontaminated 1.5 ml eppendorf tube. To homogenize add 1ml TRIZOL reagent and 10  $\mu\text{l}$  of a 20 mg/ml concentration of glycogen\* and homogenize.
  3. Add cell lysate to the tubes containing the pre-spun Phase lock tubes and incubate for 5-10 minutes @  $15 - 30^{\circ}\text{C}$  to allow complete dissociation of nucleoprotein complexes.
  4. **PHASE SEPARATION:** Add 0.2ml of chloroform-IAA per 1 ml of TRIZOL Reagent. Shake vigorously by hand for 15 seconds and incubate them @  $15-30^{\circ}\text{C}$  for 2-3 minutes.
5. Centrifuge the samples at no more than 12,000 x g for 15 minutes @  $4^{\circ}\text{C}$ . After centrifugation the mixture should separate into a lower red , phenol chloroform phase, an intermediate phase, and a colorless upper aqueous phase (RNA will be in the upper aqueous phase). The volume is typically about 60% of the aqueous phase.
6. **RNA PRECIPITATION:** Transfer the aqueous phase to a fresh tube, and save the organic phase if isolation of protein or DNA is desired. Precipitate the aqueous phase by mixing with isopropyl alcohol. Use 0.5 ml of isopropyl per 1 ml of TRIZOL reagent. Mix samples by repeated inversion. Allow samples to incubate at  $-20^{\circ}\text{C}$  for 10 minutes, then spin at 12,000 x g for 10 minutes @  $4^{\circ}\text{C}$ .
7. Decant the supernatant by dumping the liquid into a wastebasket (don't worry the RNA will remain attached to the inside of the tube).
8. **RNA WASH:** Add 1 ml of 75% EtOH per 1 ml TRIZOL Reagent used. Mix samples to dislodge the pellet and wash salts.
9. Centrifuge samples at 7,500 x g for 5 minutes @  $4^{\circ}\text{C}$ .
10. Carefully decant the alcohol and remove all the residual alcohol. Dissolve the pellet in 100 $\lambda$  RNase free ddH<sub>2</sub>O. Incubation @  $55-60^{\circ}\text{C}$  for 10 minutes may be helpful in dissolving this solution. RNA should be stored at  $-80^{\circ}\text{C}$  after this step.

#### Necessary Reagents

2ml phase lock tubes  
1.5 ml eppendorf tubes  
1.5 ml teflon tissue homogenizer  
TRIZOL RNA extraction reagent  
20mg/ml glycogen \*optional  
Choroform/isoamyl alcohol (24:1)  
Isopropyl Alcohol  
75% EtOH  
Rnase free water

---

SUPERSCRIPT II  
5X First-Strand Buffer

Note\* typical yields that I have observed are ~ 1µg/µl or 100 µg total for two whole rat eyes.

Sample calculation:

Concentration = OD260 x dilution factor x 40µg/ul

Concentration = .226 x 300/3 x 40 = 904µg/ml or 0.9µg/µl

11. Obtain the OD<sup>260/280</sup> readings by spectrophotometry to assess the purity and concentration of this extract.

12. **DNaseI Treatment for cDNA synthesis:** Mix together

- 1µg RNA sample
- 1 µl 10x Dnase Rxn Buffer
- 1µl DnaseI, amp grade, 1U/µl
- DEPC-treated H<sub>2</sub>O to 10µl

12. Incubate for 15 minutes @ RT

13. Inactivate DNaseI with 1µl 25mM EDTA solution @ 65°C for 10 min. Store RNA at -80°C until needed.