

SponSorS

The Gibbs Conference would like to thank the following companies for their support of the 26th Annual meeting.

Aviv Biomedical

Chreston

Jasco, Inc.

Micromath

GE Healthcare/Microcal

PTI


Olis, Inc.

Horiba Scientific

ISS, Incorporated

University of Alabama at Birmingham, Department of Chemistry





MicroCal™ ITC systems

Affinity – Specificity – Stoichiometry

Deeper insights into binding mechanism for confident conclusions

- Direct and automated biochemical technique
- Enable increased understanding of low and high affinity interactions
- Provide deeper insight of biological activity for a wide range of research applications



Make confident decisions sooner, visit
www.gelifesciences.com/microcal



GE, imagination at work, and GE monogram are trademarks of General Electric Company. MicroCal is a trademark of GE Healthcare companies. © 2011 General Electric Company – All rights reserved. First published Aug. 2011 GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden

HORIBA

Scientific

Sometimes combining two things gives you something much better.

Now you can have the world's fastest fluorometer and get a UV-VIS free!

HORIBA Scientific combined the world's only CCD-based benchtop spectrofluorometer with a UV-VIS spectrophotometer and created the Dual-FL. The Dual-FL features spectral rates as fast as 80,000 nm/s, and a water Raman S/N sensitivity greater than 20,000:1 RMS. Collect incredibly rapid kinetics data, or complete excitation emission matrix (EEMs) in seconds.



Oh yeah. It also has a built-in UV-VIS spectrophotometer. Use it as your lab UV-VIS, or collect simultaneous dual-mode fluorescence and absorbance data to correct inner-filter effects in your spectra.

Contact your nearest HORIBA representative to see just how fast the Dual-FL really is, and get your UV-VIS data for free.

dualfluor.com
email: adsci-tshow@horiba.com

Solutions for:

Elemental Analysis	Forensics
Fluorescence	Optical Components
Ellipsometry	Particle Characterization
Raman	



chronos



Microscopy

- FLIM, FRET
- RICS
- FCS, PCH
- Polarization

Fluorescence

- Time-Resolved
- Steady-State
- TCSPC
- Phase & Modulation

Chronos Features:

- TCSPC and/or Phase & Modulation
- Dual emission channels for T-format polarization and kinetics
- Fully computer-controlled
- Programmable sequence of measurements
- Upgradable to steady-state

ChronosBH—ISS technology combined with Becker & Hickl electronics for precise TCSPC measurements using laser diodes and LEDs.



Learn more at www.iss.com

When you Think CD... Think **JASCO**


The World's Leading Circular Dichroism Spectrometer, the J-815

JASCO's unparalleled optical performance and optionally available measurement modes make the J-815 a true "Chiro-optical Spectroscopy workbench." With the ability to measure CD, Absorbance, and Fluorescence versus wavelength and temperature on up to six samples at a time, the J-815 can provide considerable information to characterize samples faster than ever before. Also available are options for stopped-flow kinetics, auto titration, high-throughput CD and many others all controlled under the newly released Spectra Manager II™ software platform.



J-815 Circular Dichroism Spectrometer

Thermal Ramping



Multiple Position Peltier Cell Holder



MPTC-490

JASCO www.jascoinc.com

JASCO, USA • 28600 Mary's Court Easton, Maryland 21601 • 800-333-5272 • sales@jascoinc.com

The **UAB Department of Chemistry** and **UAB Center for Structural Biology** are pleased to support the 26th Annual Gibbs Conference on Biothermodynamics

UAB Department of Chemistry

College of Arts & Sciences

Graduate Studies in the Department of Chemistry

- Drug discovery, medicinal chemistry, organic synthesis
- Biophysical chemistry, structural biology
- Inorganic and organometallic chemistry
- Polymer synthesis, characterization, assembly
- Advanced materials design and development
- Computational chemistry
- Chemical education



Contact Information: Department of Chemistry Ph.D. Program
University of Alabama at Birmingham
901 14th Avenue South, Birmingham, AL 35294-1240
voice: (205) 934-4747, email: knightsen@uab.edu
Visit us on the web at www.uab.edu/chemistry



Interdisciplinary Graduate Program in Structural Biology

Graduate students in UAB's world class structural biology program experience a broad array of biophysical techniques used to address challenging biological problems at the atomic/molecular level. Our structural biology faculty are drawn from departments across the entire campus and are a truly interdisciplinary group.

Students enrolled in the structural biology program can participate in a range of structural biology courses and gain expertise on state-of-the-art equipment including high brilliance in-house X-ray generators, dedicated access to two powerful X-ray beamlines at the Argonne National lab's Synchrotron facility, high-throughput nanocrystallization systems, ultra-high resolution NMR instrumentation, large computer clusters that support macromolecular modeling and in silico drug discovery, electron microscopes for high resolution cryo-electron microscopy and tomography, and each of the major mass spectrometry technologies to support most types of structural biology-related projects.



Contact: Ms. Dianne Vickers
Program Coordinator, UAB
Biochemistry & Structural
Biology Theme
720 20th Street South, Kaul 540
Birmingham, AL 35294-0024
Tel: 205-934-0580
Email: dvickers@uab.edu